

## DIAZOTROPHS\_CO2 dinitrogen fixation and *nifH* analytical methods.

Rates of  $^{15}\text{N}_2$  fixation were assessed using the  $^{15}\text{N}$  isotopic tracer technique. Whole seawater samples from six discrete depths (5, 25, 45, 75, 100, and 125 m) were withdrawn into acid-washed duplicate 4.3 L polycarbonate bottles and sealed with caps fitted with a silicone septa. Three mL of  $^{15}\text{N}_2$  gas (98 %; Isotech Laboratories ©) was injected into each polycarbonate bottle through the silicone septum using a gas-tight syringe fitted with a stainless steel needle. The bottles were attached to a free drifting array, deployed before dawn and incubated at *in situ* light and temperature for ~24 h. After recovery of the array, the entire volume of each bottle was filtered onto pre-combusted glass microfiber (GF/F) filters (25 mm diameter) for whole seawater  $\text{N}_2$ -fixation estimates. All GF/F filters were placed on combusted foil pieces in Petri dishes and stored frozen at  $-20\text{ }^\circ\text{C}$  for transport. At the shore based laboratory, GF/F filters were dried for 24 h at  $60^\circ\text{C}$ , pelleted and  $^{15}\text{N}$  assimilation on each filter was analyzed with an elemental analyzer-isotope ratio mass spectrometer (Carlo-Ebra EA NC2500).

For the determination of cyanobacterial *nifH* gene dynamics, whole seawater samples were collected from the same six discrete depths where  $\text{N}_2$  fixation was measured into clean, acid-washed 4 L polycarbonate bottles and ~2 L was filtered through 25 mm diameter 0.22  $\mu\text{m}$  Durapore® filters (Millipore) for DNA. These filters were stored in microcentrifuge tubes containing 500  $\mu\text{L}$  DNA lysis buffer (20 mM Tris-HCL, pH 8.0; 2 mM EDTA, pH 8.0; 1.2% Triton X and 20 mg  $\text{mL}^{-1}$  lysozyme), quickly frozen in liquid nitrogen and kept at  $-80^\circ\text{C}$  until extraction. DNA samples were extracted using the DNeasy kit (Quiagen®) and RNA samples using the RNeasy Mini Kit (Quiagen®) following the manufacturer's protocol. Determinations of diazotroph gene abundances (gene copies  $\text{L}^{-1}$ ) were assessed by the quantitative polymerase chain reaction (qPCR) technique. Six different *nifH* genes (*Trichodesmium* spp., the uncultivated cyanobacteria termed "Group A", *Crocospaera* sp. and three uncultivated heterocystous cyanobacteria termed "Het-1, Het-2 and Het-3") were amplified in duplicate reactions using a TaqMan 5'-fluorogenic exonuclease qPCR assay approach. Each 25  $\mu\text{L}$  TaqMan PCR reaction was carried out using an Applied Biosystems 7300 real-time thermocycler with each sample containing 12.5  $\mu\text{L}$  TaqMan 2x Universal PCR MasterMix (Applied

Biosystems), 7.25  $\mu\text{L}$  nuclease-free, molecular-biology grade water (Sigma), 1.0  $\mu\text{L}$  each of 10  $\mu\text{M}$  forward and reverse phylotype-specific *nifH* primers (Church et al. 2005a,b), 0.5  $\mu\text{L}$  of 10  $\mu\text{M}$  fluorogenic probe, 0.75  $\mu\text{L}$  of 100x Bovine Serum Albumin (New England Biolabs), and 2  $\mu\text{L}$  of DNA extracted with the methods described previously. A six point (ranging  $10^4$  to  $10^9$  genes per reaction) standard curve was created by serial dilutions of plasmids containing the *nifH* gene inserts with a 10-fold dilution.